


# Site Planning Guide

## Mobile CT Systems

### 48' L x 8'-6" W x 13'-6" H USA Unit



### 42' L x 8'-6" W x 13'-6" H USA Unit



© 2020 Advanced Mobility Specialty Vehicles, All Rights Reserved

This manual contains confidential information of Advanced Mobility Specialty Vehicles.

This manual may be used only by you, and only for the purposes for which it was intended. You may not disclose this manual or the confidential information it contains outside of your company.

If you wish to copy any part of this manual, or to use it other than as described above, you must contact Advanced Mobility Specialty Vehicles seeking permission to do so.

You may not copy it or any part of it without the written permission of Advanced Mobility Specialty Vehicles.

As part of Advanced Mobility Specialty Vehicles' on-going program to improve its products and service, Advanced Mobility Specialty Vehicles reserves the right to implement product changes and disseminate changes in design and service information without notice or recourse.

For questions regarding the Operation or Service of this unit, please call Advanced Mobility Specialty Vehicles at 1-708-235-2800.

## List of Revisions

### Revisions

<b>Date</b>	<b>Revision Number</b>	<b>Revised by</b>	<b>Description</b>
<b>3/17/2020</b>	<b>01</b>	<b>Ramon Lopez III</b>	<b>Generalized document for all CT systems, updated information</b>

### Notice

In accordance with our policy of product development, Advanced Mobility Specialty Vehicles reserves the right to make changes in the equipment, design, specifications, and materials of the product described herein. If there are any inconsistencies between this manual and the mobile unit that inhibit serviceability, please contact Advanced Mobility Specialty Vehicles for assistance.

This manual is provided in the mobile unit. The documentation package should be kept in the mobile unit at all times.

Any problems or questions related to the components or systems covered in this manual please direct to:

#### **Advanced Mobility Specialty Vehicles**

500 Crossing Drive, Unit A  
University Park, Illinois  
60484  
USA  
(001) 708-235-2800  
(001) 708-235-2002 FAX

<http://www.amstcorp.com>

# Table of Contents

<b><u>Introduction.....</u></b>	<b>5</b>
Warnings & Safety Alert Conventions .....	5
<b><u>Support and Service Pad Requirements .....</u></b>	<b>6</b>
Trailer Weight .....	6
Recommended Support and Service Pad Requirements .....	6
Pad Depth.....	7
Pad Levelness .....	7
Electro Magnetic Interference .....	7
Vehicle Access .....	7
Recommended Attachment to the Facility.....	7
Swing Clearance Note.....	7
<b><u>Radiation Shielding Requirements .....</u></b>	<b>9</b>
Radiation Shielding.....	9
Radiation Field Information .....	9
<b><u>Customer Power Requirements .....</u></b>	<b>9</b>
Lockout/Tagout.....	9
Electrical Service .....	9
Configuration .....	9
Load Regulation at Line Frequency .....	10
Frequency.....	10
Phase Balance .....	10
Maximum Voltage Variation .....	10
Connector Type .....	10
Customer Facility.....	10
Input Power .....	10
Power Source Monitoring (Facility Only) .....	10
<b><u>Mobile Grounding Requirements .....</u></b>	<b>11</b>
Special Ground Note .....	11
<b><u>Telephone and Data Service Requirements .....</u></b>	<b>13</b>
Telephone Service.....	13
Data Service .....	13

**Water Requirements..... 13**

Water Supply Tank .....13

Portable Sink (optional) .....13

**Table of Figures**

Figure 1: Ground Connection.....11

Figure 2: Plan Layout .....15

Figure 3: Right Side Elevation .....16

Figure 4: Left Side Elevation .....18

Figure 5: Radiation Shielding Plan View .....20

Figure 6: Stair / Lift / Wall Elevation .....22

Figure 7: Russellstoll Receptacle, Service Disconnect.....26

Figure 8: Russellstoll Receptacle Chart.....27

Figure 9: Turning Requirements .....29

## Introduction

The purpose of this document is to provide the basic information needed for site planning. For specific information not contained in this document, please contact Advanced Mobility Specialty Vehicles.

The mobile unit requires sufficient room to be maneuvered and positioned for setup and takedown. The mobile unit has many storage compartments and service doors that require access during these procedures as well as during operation. The expanding wall sections, patient lift, entry stairs, and optional platform require additional space on the right side of the mobile unit. Refer to the drawings provided for actual locations of doors, patient lift, and stair sizes and locations.

## Warnings & Safety Alert Conventions

Three types of statements are used throughout this document to warn the operator of potential situations. Always read these statements carefully and take the appropriate safety precautions to ensure a safe environment for all personnel and all property. The statements are as follows:



This type of notice indicates a potentially hazardous situation, which if not avoided, could result in injury or death to the operator of the mobile unit.



This type of notice indicates a potentially hazardous situation, which if not avoided, could result in irreparable damage to the mobile unit.



This type of notice is meant to inform the operator of useful information.

## Support and Service Pad Requirements

### IMPORTANT

If other modalities utilize the same support pad, it is recommended that non-ferrous reinforcement materials be used for pad reinforcement.

### IMPORTANT

OEM of the scanning equipment must approve plans for pad construction.

The following is a list of recommendations and requirements for a concrete support and service pad. However, due to varying site conditions, the actual pad design should be prepared by an appropriately licensed structural or architectural engineer.

### Trailer Weight


The weight of the trailer should be considered in the design of the support and service pads. The overall weight of the 48' trailer is approximately 48,520lbs. The weight on the rear axles is approximately 32,780 lbs. The weight on the king pin is approximately 15,740 lbs.

Note: Weight differs for 42' trailer.


### Recommended Support and Service Pad Requirements

The measurements for the recommended support and service pad are as follows, 23' x 57'-5/8" for the 48' trailer and 20'-8" x 50'-3/4" for the 42' trailer.. The cross hatching as shown below represents the recommended support pad.

(42')



and



(42')

### **Pad Depth**

Recommendations for the width and length of the pad are given above. Based upon the existing site conditions, the depth should be determined by a local contractor.

### **Pad Levelness**

In order to ensure proper operation of the CT systems, the support pad(s) must be level and the set-up deviation must not exceed .125" in 10'-0.

### **Electro Magnetic Interference**

The ambient static magnetic field within the region of the gantry should not exceed 1 Gauss ( $10^{-4}$  Tesla) peak at the detector.

### **Vehicle Access**

A firm, level surface is required around the mobile unit in order to provide access to the site, patient access to the mobile unit, and servicing of the mobile unit.

### **Recommended Attachment to the Facility**

An inflatable air bag or soft seal is recommended at the point of connection from the unit to the facility. Fixed or solid connections may hinder imaging quality. Contact Advanced Mobility Specialty Vehicles or the local OEM of the scanning equipment representative prior to construction if the proposed connection varies from the recommended.

### **Swing Clearance Note**

Please verify the actual dimensions of the rearmost projections on the cab of your tractor to the centerline of tandem suspension or centerline of the fifth wheel plate on your tractor. Refer to [Figure 9: Turning Requirements](#) for proper tractor sizing information.





**AMST**  
500 Crossings Drive, Unit A  
University Park, IL 60484  
(708) 235-2800 main number  
(708) 235-2002 fax

*Custom Specialty Solutions*  
CT-D20-01  
Page 9 of 30  
An ISO9001:2015 Certified Company

**Page intentionally left blank.**

# Radiation Shielding Requirements

## Radiation Shielding

### **IMPORTANT**

Radiation exposure limits must be in accordance with all local, state, and federal requirements. It is the responsibility of the customer to perform a proper radiation survey in order to determine the exclusion zone.

Care should be taken when determining a site location. Factors such as shielding design, proximity to buildings, and occupancy of the surrounding areas must be considered. An exclusion zone around the mobile unit may be necessary. Refer to [\(42'\)](#)

for additional information.

## Radiation Field Information

It is the responsibility of the customer to ensure a safe environment with respect to the radiation field. Due to radioactivity levels associated with pet patient handling and diagnostic procedures used in CT scanning, an exclusion zone must be maintained while in use.

Customer must contact their local Radiation Safety Operation Official for the federal, state, and local guidelines and must comply with these safety requirements.

Operator needs to make their own exposure dose measurements to include radiation from patients when determining the outside "Keep Away Zone" (chained-off area).

# Customer Power Requirements



It is the operator's responsibility to verify that the shore power receptacle is of the same type and voltage as the connection that is supplied by Advanced Mobility Specialty Vehicles. Failure to do this can result in injury or death to the operator of the mobile unit as well as irreparable damage to the mobile unit.



The standard connector for the unit is a Russellstoll DS2504MP 480V 200A Plug. If an existing site currently implements a different connector or connector configuration, please contact Advanced Mobility Specialty Vehicles in order to arrange for a compatible power connector before the unit leaves the facility.

## Lockout/Tagout

A Lockout/Tagout provision in accordance with OSHA Standard 1910.147 is required. The facility shore power disconnect device must be located within 40'- 0" of the unit and must provide for an effective lockout/tagout to facilitate safe service and maintenance of the unit.

## Electrical Service

A single electrical power source is required for operation of the CT system. 3/N/PE AC 480V service fused at 150 amperes.

## Configuration

Three phase, five wire, wye connection, with neutral and ground. (5 wire 3/N/PE AC 480V)

## **Load Regulation at Line Frequency**

Wires are to be sized such that the line voltage drops from the power source to the mobile unit is less than 2.5% of the nominal voltage for the rated load of the mobile unit.

### **Frequency**

60Hz  $\pm$ 2.0Hz.

### **Phase Balance**

The phase balance is 3% maximum of lowest phase-to-phase voltage.

### **Maximum Voltage Variation**

The maximum voltage variation is +11% / -4% from a nominal steady state (under the worst case conditions of line voltage).

### **Connector Type**

The mobile unit is supplied with a 50'-0" power cable and male conductor. Unless otherwise specified, the connector type is a Russellstoll DS2504MP000/DF2032 480V 200A rated plug.

### **Customer Facility**

The customer facility must have the matching receptacle as specified in [Figure 7: Russellstoll Receptacle, Service Disconnect](#) and [Figure 8: Russellstoll Receptacle Chart](#). Unless otherwise specified, the receptacle type to be used must be a Russellstoll DF2504FRAB0 female connector.

### **Input Power**

- Frequency: 60Hz  $\pm$ 2.0Hz
- Regulation: Load regulation must not exceed 2.5%.
- Phase Imbalance: The difference between the highest line-to-line voltage and lowest line-to-line voltage must not exceed 3% of the lowest line-to-line voltage.

### **Power Source Monitoring (Facility Only)**

NOTE: Perform a power audit first.

A power analyzer should be used to check the proposed Mobile Siemens Truepoint PET/CT Series facility site power for average line voltage, surges, sags, reclosures, impulses, frequency and microcuts. A period that includes two weekends should be used to simulate several days of normal use. Analysis of the data and site history of any previous power problems with other X-ray systems or computer installations should be reviewed with your power and ground representative. Verify "brown-out" (low voltage) conditions, which may occur during summer months, will not exceed the allowable range.

Some analyzer models that are suitable for power monitoring are:

- Dranetz Model 658
- Dranetz Model 656A
- BMI 3630
- RPM

## Mobile Grounding Requirements

### Special Ground Note

The mobile unit must have an earth driven ground rod within 5'-0" of the facility power receptacle. A grounding cable of a minimum #1/0 AWG must be connected between the grounding rod and the grounding pin of the facility power receptacle. Another cable to be kept as short as possible may also be connected between the ground stud on the Automatic Transfer Switch and an earth driven ground rod. See [Figure 1: Ground Connection](#) below. A separate grounding conductor must still be run with the phase conductors to the source of the power from the grounding pin of the hospital power receptacle in accordance with NEC 2002 Article 250-24.




[Figure 1: Ground Connection](#)

## MOBILE GROUNDING REQUIREMENTS

**NOTE:**

- ALL WORK TO BE DONE IN ACCORDANCE WITH NATIONAL AND LOCAL ELECTRICAL CODES INFORMATION SHOWN HERE IS ONLY A RECOMMENDATION. MUST BE VERIFIED FOR SITE NATIONAL AND LOCAL CODES.
- GROUND WIRES INSIDE ENCLOSURES TO BE TAPED GREEN FOR ENTIRE VISUAL LENGTH FOR IDENTIFICATION.
- ▷ MAIN BONDING JUMPER BETWEEN GROUNDED (NEUTRAL) CONDUCTOR AND EQUIPMENT GROUNDING CONDUCTOR TO BE PROVIDED IN FACILITY SERVICE EQUIPMENT AND DOWNSTREAM AT SEPARATELY DERIVED SYSTEM TRANSFORMER SECONDARY AS SHOWN.



**GROUNDING**

The ground for our system shall originate at the system power source, i.e., transformer or first access point of power into a facility, and be continuous to our system power disconnect in the room. This ground can be spliced with "High Compression Fittings" and should be terminated at each distribution panel it passes through. When it is broken for a connection to a panel, it shall be connected into an approved grounding block with the incoming and outgoing ground in this same grounding block, which is then connected to the steel panel never using the steel panel as the block or other material of the panel as the block.

The connection at the power source shall be at the grounding point of the "Neutral-Ground" if a "Wye" transformer is used, or typical grounding points of a separately derived system. In the case of an external facility, it shall be bonded to the facility ground point at the service entrance.

**GROUNDING WIRE**

The ground wire shall be copper wire with a minimum of AWG 1/0 or the same size as the power feeders whichever is larger. This means that if there is a primary feeder to a distribution panel of 500 MCM with a secondary feeder to our system of AWG 1/0 wire, the ground to the distribution panel shall be 500 MCM with and AWG 1/0 to our system. The ground wire impedance from our system disconnect, including the ground rod shall not have an impedance greater than 2 ohms to earth as measured by one of the applicable techniques described in Section 4 of ANSI/IEEE Standard 142-1982.

Harry E. Rauworth  
Debra C. Balis  
April 22, 1999

# Telephone and Data Service Requirements

## Telephone Service

The mobile unit is supplied with up to two (2) telephone connections. The connector type that is used is a Hubbell model PH-6595 (inlet) with a model PH-6624 connector body.

The customer is required to purchase and install up to two (2) Hubbell all weather telephone connections, model PH-6597 for use at the site.

Three Hubbell model PH-6599 telephone-connecting cables are included with the mobile unit. The cables measure 50'-0" in length.

## Data Service

The mobile unit is supplied with six (6) data line connections that utilize RJ-45 outlets.

The customer is required to purchase the data connection cables for use with the data line connections. The data line connections require a 50'-0" CAT-5E cable with RJ-45 connections.

# Water Requirements

## **IMPORTANT**

During winter conditions, provisions must be made to ensure that water lines do not freeze because of weather conditions.

## Water Supply Tank

A 35-gallon water supply tank is located on the left side of the 48' mobile unit in the underbody compartments, which supplies the HVAC system.

A 20-gallon water supply tank is located in an underbody storage compartment on the road side of the 42' mobile unit, which supplies the HVAC system.


The water supply tank can be filled from within the compartment by using the supplied adapter or from the exterior of the mobile unit by using the connection on the underbody compartment door and the supplied hose.

The drain for the water supply tank is located below the underbody compartment door. The drain valve is located in the underbody compartment.


## Portable Sink (optional)

An optional portable self-contained sink is available. Within the portable sink assembly is the water supply and wastewater tank.

**Figure 2: Plan Layout  
(48')**




(42')






**Figure 3: Right Side Elevation**

(48')




(42')




**Figure 4: Left Side Elevation**

(48')



**Water Requirements (continued)**




***Water Requirements (continued)***

**Figure 5: Radiation Shielding Plan View**

(48')


**Water Requirements (continued)**



**NOTE:**  
ALL WALLS IN SCAN ROOM SECTIONS TO BE 1/16" LEAD LINED ALL EXPANDING WALLS 5/64"  
LEAD IN SCAN ROOM WALLS TO BE LEAD 84" HEIGHT ABOVE FINISHED FLOOR  
ALL FLOOR, WALL AND CORNER SEAMS TO HAVE OVERLAPPING LEAD

**Water Requirements (continued)**

(42')




***Water Requirements (continued)***

**Figure 6: Stair / Lift / Wall Elevation**


**(48')**




**Water Requirements (continued)**



**REAR ELEVATION VIEW**





**FRONT ELEVATION VIEW**




(42')

**Water Requirements (continued)**




**Figure 7: Russellstoll Receptacle, Service Disconnect**



**Water Requirements (continued)**

**Figure 8: Russellstoll Receptacle Chart**



***Water Requirements (continued)***

**Figure 9: Turning Requirements**

